PELOBiotech

www.pelobiotech.com

PELOBiotech GmbH

Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0

fax + 49 (0)89 517 286 59-88

Cellovations® Endothelial Cell Growth Medium enhanced Supplement Kit

Product Description:
Cellovations® Endothelial Cell Growth Medium Kit enhanced is formulated for optimal proliferation of endothelial cells which are cultivated under conditions described in the instructions for the used endothelial cell type. The mixture has all necessary growth stimulating elements for best results. It will give your cells optimal conditions in combination with our basal medium (PB-BH-100-9806)

The supplement is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage of the supplement:
Please note that the supplement has to be stored frozen at -20°C.

Preparing and Stability:

Please rinse the vials with 70% ethanol as protection. Remove ethanol excess before opening the cap of the vials. Please be careful not to touch the vials interior threads with your fingers. Also do not warm up the supplement over 37°C. All this protects your supplement against contamination or degradation of essential components.

To cultivate large vessel endothelial cells under optimal conditions please add our supplement kit components to a basal medium; best results will be with our basal medium (Cellovations® Endotheleial Cell Basal Medium, PB-BH-100-9806)!

Product Information

- Cat. No.: PB-SH-100-2199
- Components:
  Cellovations® Endothelial Cell Growth Medium enhanced Supplement Kit
  - Glutamin
  - FCS
  - bFGF
  - hEGF
  - VEGF
  - R3-IGF-1
  - Hydrocortison
  - Heparin
- *Optional:*
  Antibiotics Solution (Gentamycin 50 µg/mL, Amphotericin 50 µg/mL)
  (Cat. No.: PB-SX-000-0099)