PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59 88

www.pelobiotech.com

PassagePack

PassagePack

Product Description:

The PassagePack contains Trypsin/EDTA (T/E), Trypsin Neutralization Solution and HBSS an optimized combination to be used to detach cells from culture vessels.The solutions are sterile filtered through a 0.2µm membrane filter.

Components:

Trypsin/EDTA Solution in HBSS, w/o Ca2+ and Mg2++ Cat# PB-LU-002-0100 (100 ml)
Trypsin Neutralization Solution in HBSS, w/o Ca2+ and Mg2+ Cat# PB-LU-002-0110 (100 ml)
Hank's Balanced Salt Solution(HBSS) (Cell culture grade) Cat# PB-LU-003-0010 (100 ml)

Procedure of use

1) Remove culture medium from culture vessels by aspiration and wash the cells with HBSS w/o Ca^{2+} and Mg^{2+} to remove all traces of serum. Remove HBSS w/o Ca^{2+} and Mg^{2+} by aspiration.
2) Dispense enough trypsin/EDTA solution into culture vessel to completely cover the monolayer of cells and place in room temperature (or 37°C incubator if necessary) for approximately 1-2 minutes or until 80% of cells are rounded up (monitored with microscope). (e.g. 6 ml per T-75 flask)
3) Remove the trypsin/EDTA solution by aspiration and return closed culture vessel to incubator for another 1 minute (hit the culture vessel firmly against bench top, all cells should be moving on the surface of culture vessels under microscope).
4) Add 5 ml trypsin Neutralization Solution to the cells as soon as possible to inhibit further tryptic activity which may damage cells.
5) Cells can be harvested by gently pipetting the cell suspension. Further dilution can be made, if required, for cell counts and/or subculturing.

Specifications

pH	7.4 ± 0.2
Cell Culture Test	Pass
Sterility	Pass

Storage: Stored at -20°C for long term storage. Once thawed, keep it in a refrigerator at 2-4°C for one month.

Quality control:

1. Sterility test (bacteria, fungi and mycoplasma)
2. pH test

Product Information

Catalogue number:	PB-900028	3 x 100 ml

THESE PRODUCTS ARE FOR RESEARCH USE ONLY

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88