

## Cellovations® Keratinocyte Growth Media

### Cellovations® Keratinocyte Growth Medium Kit enhanced, defined & ACF (ready-to-use)

- optimized growth conditions for human keratinocytes

**Product Description:**

Cellovations® Keratinocyte Growth Medium Kit enhanced, defined & ACF is a defined low calcium (0,07 mM) medium formulated for optimal growth of human keratinocytes. It contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals. This medium is animal component free (ACF). **The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.**

**Storage Conditions:**

Store medium at -20°C.

Protect from light for a longer time period.

**Preparing and Stability:**

Thaw the medium at 37°C. After thawing please rinse the bottle and the tube with 70% ethanol as protection. Remove ethanol excess before opening the cap of the bottle.

Please be careful not to touch the bottle interior threads with your fingers. Please avoid long term exposure to light because some components are light sensitive. Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components.

---

**Product Information:**

- Cat. No.: PB-C-MH-300-1099-d
- Components: 500 ml ready-to-use Medium

MSDS are available on request!

---

**ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88**