

Pericytes Growth Medium

Pericytes Growth Medium (Ready to use)

Product Description:

Pericyte growth medium is consisted of a high glucose basal medium and pericyte growth supplements (including fetal bovine serum, growth factors). Due to the nature of those elements, certain level of lot-to-lot variability is unavoidable.

Each lot of this product is tested by culturing human pericytes to ensure its full activity and sterility (This product contains antibiotics: penicillin and streptomycin).

Storage of the medium and supplements

The Pericyte growth medium should be stored at 4°C. A change in color or appearance of precipitate may indicate deterioration.

Shelf Life

The Pericyte growth medium ready to use expires after at least 2 months from the date of receipt under proper storage condition.

Preparing and Stability:

Before using, transfer just required amount (ml) of PC-Growth Medium for your experiment from the original bottle and prewarm it to 37°C in a water bath or an incubator (Avoid frequent temperature change to the entire bottle of medium).

We recommend changing medium every 2 days for regular cell culture.

Product Information

- Catalogue number: PB-MH-031-4000 Pericyte Growth Medium (Ready to use)
500 ml

THIS PRODUCT IS FOR RESEARCH USE ONLY and NOT FOR DIAGNOSTIC OR
THERAPEUTIC APPLICATIONS

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88