

Bladder Epithelial Cell Growth Medium

Bladder Epithelial Cell Growth Medium Kit, serum-free

- optimized growth conditions for human bladder epithelial cells

Product Description:

Cellovations® Bladder Epithelial Cell Growth Medium Kit is a complete medium formulated for optimal growth of bladder epithelial cells. The medium contains amino acids, minerals, vitamins and organic compounds and supplements. It does not contain antibiotics / anti-mycotics. **FCS, 50 ml, must be added.**

The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage of the medium:

Please note that the basal medium has to be stored at 4°C up to 8°C but never should be frozen. The supplements have to be stored at -20°C. Protect from light for a longer time period.

Preparation and Stability:

Thaw the supplement KIT at 37°C before adding to the basal medium. After thawing please rinse the bottle and the tubes with 70% ethanol as protection. Remove ethanol excess before opening the caps of the tubes. After adding the supplements to the basal medium slew carefully to get a homogenous medium.

Please be careful not to touch the tubes or bottle interior threads with your fingers. **Please avoid long term exposure to light because some components are light sensitive.** Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components.

Product Information

Catalogue number:	PB-MH-355-0090-SF
Components:	500 ml Basal Medium – Supplement KIT – Antibiotics Solution
(optional)	0.5 mL Insulin-Transferrin-Selenium (ITS) 0.5 mL EGF 0.5 mL Hydrocortisone 5.0 mL L-Glutamine 5.0 mL Antibiotic-Antimycotic Solution 10.0 mL Epithelial Cell Supplement 50.0 mL FBS have to be added by customer

Orders can be placed at sales@pelobiotech.com or by fax at + 49 (0)89 517 286 59 88