

# Freezing Medium, defined, animal-free

## Freezing Medium

### Product Description:

PELOBIOTECH Freezing Medium is a 2 x concentrate, cryo preservation medium containing 20% DMSO and all the necessary components to protect cells during freezing and thawing. This freezing medium is chemically defined and animal component free. It is ready to use.

Formulated for optimized preservation of cells during freezing and thawing. Provided as an 2 x concentrate, use as a 1 x freezing medium. For freezing and thawing protocols, please see our cell freezing protocol.

### Special Note:

We recommend to use our Freezing Medium together with a special freezing container like a CoolCell or Mr. Frosty.

### Storage:

Store freezing medium at 4 °C in the dark.

### Quality Control:

Growth and morphology of cells tested after thawing using the recommended protocol. Free of bacteria, fungi and mycoplasma contamination.

---

### Product Information

Catalogue number:	PB-PACRY-50
Quantity:	50 ml, 2x concentrate

**For research use only. Not for use in therapy or diagnostics.**

**ORDER UNDER [sales@pelobiotech.com](mailto:sales@pelobiotech.com) or FAX + 49 (0)89 517 286 59 88**