PELOBIOTECH GmbH
Behringstr.10
82152 Planegg Germany
info@pelobiotech.com
www.pelobiotech.com
tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59-88

PELOBiotech
www.pelobiotech.com

Human Lung Microvascular Endothelial Cells

Human Lung Microvascular Endothelial Cells
Cat#: PB-CH-150-4011

Product Description:
Human Lung Microvascular Endothelial Cells (HLMECs) (PB-CH-150-4011) are isolated from healthy normal human lung tissue from an adult donor.
The cells are shipped in cryopreserved culture (the cells are provided @ passage 3). Endothelial Growth Medium MV enhanced Kit (contains 5% serum and growth supplements, Cat# PB-MH-100-4099) is recommended for cell culture and these cells have an average additional population doubling levels >**15** when cultured following the detailed protocol.

Storage of cells:

When you receive the cells please store them immediately in liquid nitrogen.

NOT FOR CLINICAL, DIAGNOSTIC OR Therapeutic Procedures. DO NOT USE IN HUMANS. Since the material is of human origin, these cells have to be treated as potentially biologically hazardous material, even if serological tests are negative.

Characterization of the Cells

Cytoplasmic VWF / Factor VIII:	**>95% positive by immunofluorescence**
Cytoplasmic uptake of Di-I-Ac-LDL:	**>95% positive by immunofluorescence**
Cytoplasmic PECAM1:	**>95% positive by immunofluorescence**
Initial Seeding Density after thawing:	2 x 10^4 cells / cm^2
Seeding Density after passaging:	1 x 10^4 cells / cm^2

HLMECs are negative for HIV-1, HBV, HCV, bacteria, yeast, fungi and mycoplasma.

Product Information

Name of Cells:	Human Lung Microvascular Endothelial Cells	
Catalogue Number:	PB-CH-150-4011	
Product Format:	cryopreserved	
Cell Number:	**>**5 x 10^5 cells / vial	
Medium:	PB-MH-100-4099	Endothelial Growth Medium MV enhanced Kit (Basal medium + Supplements).

HLMEC´s are for research use only.

ORDER UNDER: sales@pelobiotech.com or FAX + 49 (0)89 517 286 59-0