PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59-88

www.pelobiotech.com

Human Gingival Epithelial Cells, pooled

Human Gingival Epithelial Cells, pooled
Cat#: PB-CH-341-7413

Product Description:
Human Gingival Epithelial Cells, pooled (HGEPp) (PB-CH-341- 7411) are pooled cells isolated from normal human gingiva tissue of 3 or more donors.

The cells are shipped in cryopreserved culture

(the cells are provided @ passage 1 to 3 (See CoA for detailed information).
Pelobiotech's Oral Epithelial Cell Growth Medium (contains no fetal calf serum only low BPE and chemically defined growth supplements, Cat# PB-C-MH-341-7499) is recommended for cell culture and these cells have an average additional population doubling levels >15 when cultured following the detailed protocol.

Storage of cells:

When you receive the cells please store them immediately in liquid nitrogen.

NOT FOR CLINICAL, DIAGNOSTIC OR Therapeutic Procedures. DO NOT USE IN HUMANS. Since the material is of human origen, these cells have to be treated as potentially biologically hazardous material, even if serological tests are negative.

Characteristics of the Cells

Initial Seeding Density:	1 x 10^4 cells / cm^2
Seeding Density:	4 x 10^3 cells / cm^2

HBEPC are negative for HIV-1, HBV, HCV, bacteria, yeast, fungi and mycoplasma.

Product Information

Name of Cells:	Human Gingival Epithelial Cells
Catalogue Number:	PB-CH-341-7413
Product Format:	cryopreserved
Cell Number:	>5 x 10^5 cells / vial

Medium: PB-C-MH-341-7499 Oral Epithelial Cell Growth Medium (Basal medium + BPE + Supplements).

HBEPC´s are for research use only.

ORDER UNDER: sales@pelobiotech.com or FAX + 49 (0)89 517 286 59-0