PELOBIOTECH GmbH
Behringstr.10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59-88

www.pelobiotech.com

Human Brain Astrocytes

Human Brain Astrocytes
Cat#: PB-CH-719-2911


Product Description:
HBA (PB-CH-719-2911) are initiated by elutriation from dissociated normal human brain cortical tissues. The cells are shipped in a frozen vial (the cells are provided @ passage 1). Astrocyte Growth Medium (PB-MH-719-0000) is recommended for cell culture and these cells have an average population doubling capacity > 10 when cultured following the detailed protocol described below.

Storage of cells:
When you receive the cells please store them immediately in liquid nitrogen.

NOT FOR CLINICAL, DIAGNOSTIC OR Therapeutic Procedures. DO NOT USE IN HUMANS. Since the material is of human origin, these cells have to be treated as potentially biologically hazardous material, even if serological tests are negative.

Characterization of the Cells

Cytoplasmic GFAP	> 98% positive by immunofluorescence
Seeding Density:	2 x 10^4 cells / cm^2



HBA are negative for HIV-1, HBV, HCV, and mycoplasma.

Product Information

Name of Cells:	Human Brain Astrocytes
Catalogue Number:	PB-CH-719-2911
Product Format:	cryopreserved
Cell Number:	>5 x 10^5 cells / vial

HBA are for research use only.

ORDER UNDER: sales@pelobiotech.com or FAX + 49 (0) 89 517 286 59-0