PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59 88

PELOBiotech
www.pelobiotech.com

Fetal Bovine Serum, EU Approved

Product Description:

Serum is the blood component obtained after coagulation and removal of cellular components. Besides serum proteins it contains growth factors, amino acids and hormones, etc. This comprehensive mix makes serum one of the most important supplements for supporting cell growth and proliferation in in vitro cell culture. Of special interest is fetal bovine serum (FBS). It is especially rich in growth factors and low in antibodies, which may influence the cell culture work.

Product Specification:

Product number: PB-FCS-EU-0100

Vol.: 100ml

Storage and shelf life: Store at ≤ -15°C.

Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended.

Once opened, store at +4°C and use within 4-6 weeks.

Shipping: Frozen on dry ice

To thaw: Keep overnight at 2-8°C, swirl gently to homogenize solution.

For research use only. Not for use in therapy or diagnostics.

Sterility Testing:

This product has been tested for and has been shown to be free of bacterial, fungal, and mycoplasma contamination.

Product Information

Catalogue number:	PB-FCS-EU-0100
Product:	FBS, EU Approved
Size:	100 ml