PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59 88

www.pelobiotech.com

Freezing Medium (FBS)

PB-LH-000-0000-50 Cell Culture Freezing Medium (FBS)

Product Description:

PELOBIOTECH Freezing Medium is a FBS based cryo preservation medium containing 10% DMSO and all the necessary components to protect cells during freezing and thawing. This freezing medium is ready to use. Formulated for optimized preservation of cells during freezing and thawing.

Size:	50ml
Storage:	< -20 °C for long term storage / keep at 2 - 8°C after thawed
Intended Use:	Cell Freezing Medium is used for cryopreservation of cells
Shelf life:	-20°C for 12 month
Protocols for Freezing Cells:	For optimum results, cells should be in log phase of growth. Please follow the recommended protocols.
Quality Control:	Cell Freezing Medium is performance tested on HUVECs. Additional standard evaluations are endotoxin, pH, osmolality and tests for the absence of bacterial and fungal contaminants.

Product Information

Catalogue number:	PB-LH-000-0000-50
Quantity:	50 ml

For research use only. Not for use in therapy or diagnostics.

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88