PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 771 59-88

www.pelobiotech.com

Cellovations® Thawing Medium for NPCs

Cellovations® Thawing Medium for NPCs (TMNPC)
- optimized thawing conditions for NPCs

Product Description:

Cellovations® Thawing Medium for NPCs (TMNPC) is a complete medium formulated for optimal thawing of NPCs.

The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage of the medium (TMNPC)

Please note that the medium has to be stored at 4°C but never should be frozen. Protect from light for a longer time period.

Preparing and Stability:

The medium should be used onced opened within two weeks.
Please be careful not to touch the tube or bottle interior threads with your fingers. Please avoid long term exposure to light because some components are light sensitive. Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components.

Product Information

• Catalogue number:	PB-MH-053-0000-T50
Components:	50 ml Thawing Medium for NPCs
Storage:	Sore at 4°C

The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88