PELOBIOTECH GmbH
Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 771 59-88

PELOBiotech
www.pelobiotech.com

Cellovations® EASY Endothelial Cell Growth Media

Cellovations® EASY Endothelial Cell Growth Medium

- optimized growth conditions for Endothelial Cells
especially larger vessels

Product Description:
Cellovations® Endothelial Cell Growth Medium EASY is a complete medium formulated for optimal growth of endothelial cells. It contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals plus 2% FCS together with the mixture of supplements.
The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage of the medium and supplements
Please note that the basal medium has to be stored at 4°C up to 8°C but never should be frozen.
The supplement has to be stored at -20°C. Protect from light for a longer time period.

Preparing and Stability:
Thaw the supplement at 37°C before adding to the basal medium. After thawing please rinse the bottle and the tube with 70% ethanol as protection. Remove ethanol excess before opening the cap of the tube. After adding the supplement to the basal medium slew carefully to get a homogenous medium.

Please be careful not to touch the tube or bottle interior threads with your fingers. Please avoid long term exposure to light because some components are light sensitive. Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components.

Product Information

- Catalogue number: PB-MH-100-2100 Cellovations® Endothelial Cell Growth Medium EASY
- Components: 500 ml Basal Medium – Mixture of Supplements – Antibiotics Solution
- Mixture of Supplements contains: 2% FCS (fetal calf serum),
EndoCGS/Heparin,
bFGF, EGF,
Hydrocortisone,
Glutamin

optional: 50.00 ng / ml Amphotericin B; 50.00 µg / ml Gentamicin

MSDS are available on request!

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88