PELOBiotech GmbH

Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59-88

www.pelobiotech.com

Cellovations® Oral Epithelial Cell Growth Media

Cellovations® Oral and Gingival Epithelial Cell Growth Medium Kit, low BPE

- optimized growth conditions for human oral and gingival epithelial cells

Product Description:
PELOBiotech's Cellovations® Oral Epithelial Cell Growth Medium Kit low BPE, is a serum free medium formulated for optimal growth of human gingival and oral epithelial cells. It contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, BPE, growth factors and trace minerals.
The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage Conditions:

- Basal Medium, at 4°C up to 8°C (never should be frozen)
- Supplements I - III and BPE at -20°C.

Protect from light for a longer time period.

Preparing and Stability:
Thaw the supplements at 37°C before adding to the basal medium. After thawing please rinse the bottle and the tube with 70% ethanol as protection. Remove ethanol excess before opening the cap of the tube. After adding the supplements to the basal medium slew carefully to get a homogenous medium.

Please be careful not to touch the tube or bottle interior threads with your fingers. Please avoid long term exposure to light because some components are light sensitive. Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components.

Product Information:

• Cat. No.:	PB-MH-341-7499
• Components:	500 ml Basal Medium – Supplement Kit
Supplement Kit contains:	Supplement I
	Supplement II
	Supplement III
	BPE

MSDS are available on request!

ORDER UNDER sales@pelobiotech.com or FAX + 49 (0)89 517 286 59 88