PELOBiotech
www.pelobiotech.com

PELOBiotech GmbH

Behringstr. 10
82152 Planegg Germany

info@pelobiotech.com
www.pelobiotech.com

tel + 49 (0)89 517 286 59-0
fax + 49 (0)89 517 286 59-88

Cellovations® Tenocyte Growth Media

Tenocyte Growth Medium classic; ready-to-use

- optimized medium to culture human tenocytes

Product Description:
Cellovations® Tenocyte Growth Medium classic is a complete medium formulated for optimal growth of tenocyte isolated out of patella or Achilles tendon. The medium contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals plus FCS together with the supplement mixture.
Our medium is consisted of a basal medium and already added tenocyte growth supplements (including fetal bovine serum and growth factors) in the ready to use version or separately in the kit version. Due to the nature of those elements, certain level of lot-to-lot variability is unavoidable

The medium is for research use only and not approved for human or animal use or for application in any in vitro diagnostic procedure.

Storage Conditions:

- Basal Medium at 4°C up to 8°C (never should be frozen)
- Supplements at -20°C.
- Ready To Use - Medium at 4°C up to 8°C (never should be frozen)

Protect from light for a longer time period.

Preparing and Stability:
Thaw the supplement kit at 37°C before adding to the basal medium. After thawing please rinse the bottle and the tubes with 70% ethanol as protection. Remove ethanol excess before opening the caps of the tubes. After adding the supplements to the basal medium slew carefully to get a homogenous medium.

Please be careful not to touch the tube or bottle interior threads with your fingers. Please avoid long term exposure to light because some components are light sensitive. Also do not warm up the medium over 37°C. All this protects your medium against contamination or degradation of essential components. Stored in the dark at 4°C the reconstituted medium will be stable for 5 weeks.

Product Information:

- Cat. No.: PB-MH-479-0000-rtu
- Components: 500 ml complete Medium ready to use
  500 ml Basal Medium – 50 ml Supplements

MSDS are available upon request.

Orders can be placed at sales@pelobiotech.com or by fax at + 49 (0)89 517 286 59 88